



Product Code . DS-GF-10364

Flasks-E Amber, Without Side Cut for Dissolution Apparatus

Description

Flasks-E Amber, Without Side Cut for Dissolution Apparatus

Specification:-

- o Material: Borosilicate Glass
- o Graduated in interval of 50 ml.
- o Smooth hemispherical base.
- o Capacity: 1000 mL
- o Diameter: 131mm
- o Height: 173mm
- o Heat Resistant
- o Chemical Resistant
- o Used in Chemical Laboratory

