



Product Code . DS-GF-10365

Flasks-L Amber, Without Side Cut for Dissolution Apparatus

Description

Flasks-L Amber, Without Side Cut for Dissolution Apparatus

Specification:-

- o Material: Borosilicate Glass
- o Color: Amber
- o Graduated in interval of 50 ml.
- o Smooth hemispherical base.
- o Heat Resistant
- o Chemical Resistant
- o Used in Chemical Laboratory

Didac Scientific Pvt. Ltd,
Martinfield Business Centre, 108 Martinfield, Welwyn Garden City , United Kingdom
Direct Contact Details ✉ sales@didacscientific.co.uk
💻 www.didacscientific.co.uk/