




Product Code . DS-GF-10374

Volumetric Flask, Pyrex

Description

Volumetric Flask, Pyrex

Specification:-

- o Material: Pyrex Glass
- o Capacity: 10 ml, 50 ml, 100 ml, 250 ml, 500 ml, 1000 ml, 2000 ml, 5000 ml
- o With Plastic Stopper
- o Graduated
- o Heat Resistant
- o Chemical Resistant
- o Color: Clear
- o Used in Chemical Laboratory

Didac Scientific Pvt. Ltd,
Martinfield Business Centre, 108 Martinfield, Welwyn Garden City , United Kingdom
Direct Contact Details ✉ sales@didacscientific.co.uk
💻 www.didacscientific.co.uk/