



Product Code . DS-GF-10418

Volumetric Flask Brush

Description

Volumetric Flask Brush

Descriptions:

- o Made of White Natural Bristles
- o Mounted on a Durable, plastic handle.
- o Flexible, curved brush adapts to the shape of flasks.
- o Diameter: 1 3/8"
- o Length: 3 1/2"
- o Used for Cleaning Volumetric Flask

Didac Scientific Pvt. Ltd,
Martinfield Business Centre, 108 Martinfield, Welwyn Garden City , United Kingdom
Direct Contact Details ✉ sales@didacscientific.co.uk
💻 www.didacscientific.co.uk/