



Product Code . DS-GB-10429

Micro Burettes, Class-A

Description

Micro Burettes, Class-A

Description:

These Class-A micro burettes are designed for precise analytical methods in micro chemistry. Tips are tapered to assure smooth flow. Manufactured from precision bore tubing assures accuracy between any two points throughout the body. Funnel top facilitates addition of reagents by pouring, or by suction from below through a rubber-stoppered tube.

- o Class A Burette made of high quality, heavy duty borosilicate glass.
- o With PTFE Key stopcock
- o Capacity: 5 ml, 10 ml
- o Blue graduations
- o Heat Resistant
- o Chemical Resistant
- o Used in Chemical Laboratory

