

Product Code . DS-GB-10431

Burette, Class-B, with PTFE Key Stopcock





Description

Burette, Class-B, with PTFE Key Stopcock

Description:

- o Class B Burette made of high quality, heavy duty borosilicate glass.
- o With PTFE Key Stopcock
- o Capacity: 5 ml, 10 ml, 25 ml, 50 ml, 100 ml
- o Graduated
- o Heat Resistant
- o Chemical Resistant
- o Used in Chemical Laboratory

Didac Scientific Pvt. Ltd,
Martinfield Business Centre, 108 Martinfield, Welwyn Garden City , United Kingdom
Direct Contact Details  sales@didacscientific.co.uk
 www.didacscientific.co.uk/