



**Product Code . DS-GB-10433**

## **Burettes with PTFE Needle Valve Stopcock**

### **Description**

---

#### **Burettes with PTFE Needle Valve Stopcock**

##### **Description:**

o Borosilicate glass body with a Rotaflow type PTFE Needle valve stopcock which ensures a leakproof seal, finer flow control, elimination of lubrication and absence of plug seizure.

o Capacity: 5 ml, 10 ml, 25 ml, 50 ml, 100 ml

o Graduated

o Heat Resistant

o Chemical Resistant

o Used in Chemical Laboratory

---

Didac Scientific Pvt. Ltd,  
Martinfield Business Centre, 108 Martinfield, Welwyn Garden City , United Kingdom  
Direct Contact Details ✉ [sales@didacscientific.co.uk](mailto:sales@didacscientific.co.uk)  
💻 [www.didacscientific.co.uk/](http://www.didacscientific.co.uk/)