



Product Code . DS-GB-10435

Burettes with detachable Stopcocks Schellbach Stripe

Description

Burettes with detachable Stopcocks Schellbach Stripe

Description:

- o Schellbach stripe for accurate menscus reading
- o Burettes made of high quality, heavy duty borosilicate glass.
- o With PTFE Key Stopcock
- o Capacity: 5 ml, 10 ml, 25 ml, 50 ml, 100 ml
- o Graduated
- o Heat Resistant
- o Chemical Resistant
- o Used in Chemical Laboratory

Didac Scientific Pvt. Ltd,
Martinfield Business Centre, 108 Martinfield, Welwyn Garden City , United Kingdom
Direct Contact Details ✉ sales@didacscientific.co.uk
💻 www.didacscientific.co.uk/