



Product Code . DS-GB-10438

Aspirator Bottles

Description

Aspirator Bottles

Description:-

- o Made of polypropene and are thus abundant lighter than glass.
- o These are autoclavable, sturdy & dependable.
- o These bottles are used for storing or dispensing water laboratory solutions.
- o The leakproof spigot screwed to the bottle body needs solely a communicate open or shut.
- o Threads ought to be completely disengage before autoclaving.
- o Capacity: 0.5 Ltr, 1.5 Ltr, 5 Ltr, 10 Ltr, 20 Ltr

Didac Scientific Pvt. Ltd,
Martinfield Business Centre, 108 Martinfield, Welwyn Garden City , United Kingdom
Direct Contact Details ✉ sales@didacscientific.co.uk
💻 www.didacscientific.co.uk/