



Product Code . DS-GB-10439

Reagent Bottles

Description

Reagent Bottles

Description:-

- o Specially designed PE Pouring rings for drip-free operations.
- o Clear scale Graduation & large marking area.
- o These complete bottles are autoclavable & can be sterilized.
- o Autoclavable PE screw cap and pouring ring.
- o These bottles are highly mechanically durable and chemical resistance.
- o Material: Borosilicate Glass
- o Capacity:
 - 10ml
 - 25ml
 - 50ml
 - 100ml

250ml

500ml

1000ml

2000ml

5000ml

10000ml



20000ml

o Heat Resistant

o Chemical Resistant

o Color: Clear

o Used in Chemical Laboratory

Didac Scientific Pvt. Ltd,
Martinfield Business Centre, 108 Martinfield, Welwyn Garden City , United Kingdom
Direct Contact Details  sales@didacscientific.co.uk
 www.didacscientific.co.uk/