



Product Code . DS-A-10451

Automatic Autoclave Swiftlock Ecofill

Description



Automatic Autoclave Swiftlock Ecofill

Description:

- Eco-friendly design – low water and low power consumption, the chamber is designed so the minimum amount of water is needed per cycle
- Variable temperature range, 100 to 138°C (0.2 to 2.4bar)
- Fully automatic once filled with water
- Stainless steel vessel
- Base tray / shelf
- Thermocouple validation gland with easy access panel
- Gauge test port
- Safety valve test cycle
- Over temperature and over pressure cut-out
- Cool interlock (80°C) - fluid cycles only
- One of the largest easy to use icon driven colour touchscreen's on any benchtop/compact sterilizer

Safety Features:

- Pressure activated door lock
- Overheat protection
- Sterilisation failure, water low and door unlocked alarms
- Two thermal cooling locks in compliance
- Audio-visual end of cycle, cycle fault, cycle interruption.

Didac Scientific Pvt. Ltd,
Martinfield Business Centre, 108 Martinfield, Welwyn Garden City , United Kingdom
Direct Contact Details  sales@didacscientific.co.uk
 www.didacscientific.co.uk/