



**Product Code . DS-A-10454**

## **Automatic Autoclave Swiftlock Front Load**

### **Description**

---

#### **Automatic Autoclave Swiftlock Front Load**

##### **Description:**

- Fully automatic once filled with water
- Grade 316L stainless steel vessel
- Variable temperature range, 100 to 138°C (0.2 to 2.4bar)

##### **Safety Features:**



- Pressure activated door lock
- Two thermal cooling locks in compliance.
- Insulated lid cover
- Over-pressurisation safety valve
- Audio-visual end of cycle, cycle fault, cycle interruption alarms
- Sterilisation failure, water low and door unlocked alarms

---

- Swiftlock' door closure, manual filling and base shelf.

- For 230V 50/60Hz single phase 13A supply.

---

Didac Scientific Pvt. Ltd,  
Martinfield Business Centre, 108 Martinfield, Welwyn Garden City , United Kingdom  
Direct Contact Details  [sales@didacscientific.co.uk](mailto:sales@didacscientific.co.uk)  
 [www.didacscientific.co.uk/](http://www.didacscientific.co.uk/)