



Product Code . DS-I-10455

Mini Incubator

Description

Mini Incubator

Description:

- Standard profile with controls at base
- Natural convection
- Solid door
- Aluminium coated mild steel interior
- Fluctuation $\pm 0.25^{\circ}\text{C}$ at 37°C
- Top vent with clip for thermometer
- Fitted with fixed shelf runners and removable chrome-plated wire grid shelves
- Steel outer case with stoved powder paint finish
- Hydraulic thermostat control and separate adjustable safety thermostat, both with locking devices and indicator lamps
- For 220-240V 50Hz single phase supplies.

Martinfield Business Centre, 108 Martinfield, Welwyn Garden City , United Kingdom

Direct Contact Details  sales@didacscientific.co.uk

 www.didacscientific.co.uk/