



Product Code . DS-I-10457

Incubator Mini 30L, INT. 70°C Metal Door 220-240V 50HZ

Description

Incubator Mini 30L, INT. 70°C Metal Door 220-240V 50HZ

Description:

- Aluminium coated mild steel interior
- Solid door
- Steel outer case with stoved powder paint finish
- Fluctuation $\pm 0.25^{\circ}\text{C}$ at 37°C
- Top vent with clip for thermometer
- Standard profile with controls at base
- Natural convection
- Hydraulic thermostat control and separate adjustable safety thermostat, both with locking devices and indicator lamps
- Fitted with fixed shelf runners and removable chrome-plated wire grid shelves
- For 220-240V 50Hz single phase supplies.

Specifications:

- Capacity: 30 litres
- Power rating: 1500 W
- Shelves supplied: 2
- Shelf positions: 2
- Weight: 20kg
- External H x W x D: 500 x 480 x 490mm
- Internal H x W x D: 240 x 360 x 350mm

Didac Scientific Pvt. Ltd,
Martinfield Business Centre, 108 Martinfield, Welwyn Garden City , United Kingdom
Direct Contact Details  sales@didacscientific.co.uk
 www.didacscientific.co.uk/