

Product Code . DS-I-10458



Incubator 70°C 50L, Solid Door 220-240V 50HZ

Description



Incubator 70°C 50L, Solid Door 220-240V 50HZ

Description:

- Wide bodied profile with controls to the side
- Fluctuation $\pm 0.25^{\circ}\text{C}$ at 37°C
- Natural convection
- Top vent with clip for thermometer
- Fitted with fixed shelf runners and removable chrome-plated wire grid shelves
- Steel outer case with stoved powder paint finish
- Aluminium coated mild steel interior
- Solid door
- Touchscreen digital control incorporating an integral overheat system with two Pt100 sensors.
- For 220-240V 50Hz single phase supplies.

Specifications:

- Capacity: 50 litres
- Power rating: 750 W
- Shelves supplied: 2
- Shelf positions: 3
- Weight: 31 kg
- External H x W x D: 580 x 730 x 480mm
- Internal H x W x D: 420 x 360 x 350mm

Didac Scientific Pvt. Ltd,
Martinfield Business Centre, 108 Martinfield, Welwyn Garden City , United Kingdom
Direct Contact Details  sales@didacscientific.co.uk
 www.didacscientific.co.uk/