



Product Code . DS-I-10459

Mini Incubator 50 Litres

Description



Mini Incubator 50 Litres

Description:

- Standard profile with controls at base
- Stainless steel interior
- Fluctuation $\pm 0.25^{\circ}\text{C}$ at 37°C
- Top vent with clip for thermometer
- Solid door
- Natural convection
- Fitted with fixed shelf runners and removable chrome-plated wire grid shelves
- Steel outer case with stoved powder paint finish
- Hydraulic thermostat control and separate adjustable safety thermostat, both with locking devices and indicator lamps
- For 220-240V 50Hz single phase.

Specifications:

- Capacity: 50 litres
- Power rating: 500 W
- Shelves supplied: 2
- Shelf positions: 3
- Weight: 26 kg
- External H x W x D: 590 x 610 x 470mm
- Internal H x W x D: 330 x 490 x 330mm

Didac Scientific Pvt. Ltd,
Martinfield Business Centre, 108 Martinfield, Welwyn Garden City , United Kingdom
Direct Contact Details  sales@didacscientific.co.uk
 www.didacscientific.co.uk/