



Product Code . DS-I-10460

Memmert Incubator

Description



Memmert Incubator

Description:

- Stainless steel interior and exterior sides
- Temperature range ambient +5°C to +80°C
- Natural convection
- With Single Display control Cockpit controller
- 10 year datalogger
- Setting accuracy: $\pm 0.1^\circ\text{C}$
- Supplied with stainless steel wire grid shelf
- Inner glass door
- Heaters located in ribs on all four sides provide uniform heating
- Digital process timer, adjustable between 1 minute and 99 days, 23 hours

Specifications:

- Capacity: 32 litres
- Shelves supplied: 1
- Weight: 44 kg
- Internal H x W x D: 320 x 400 x 250mm

Didac Scientific Pvt. Ltd,
Martinfield Business Centre, 108 Martinfield, Welwyn Garden City , United Kingdom
Direct Contact Details  sales@didacscientific.co.uk
 www.didacscientific.co.uk/