



Product Code . DS-I-10461

CO2/RH Incubator Memmert

Description

CO2/RH Incubator Memmert

Description:

- Adjustable humidity control from 80 to 93% RH, accuracy $\pm 2.5\%$
- Digital CO2 control from 0 to 20%
- Stainless steel interior and exterior
- Automatic start-up
- Sterilisation programme: 60 minutes at 180°C
- Classified as a Class IIa medical device for IVF work
- Adjustable temperature control from ambient + 5 to +50°C, $\pm 0.5^\circ\text{C}$ (chamber variation)
- Twin Display control Cockpit multifunctional PID controller with two high-definition TFT colour displays, timer and fault indicators, ethernet LAN and USB output interfaces

Specifications:

-
- Capacity: 107 litres
 - Temperature range: Ambient +5 to +50 deg C
 - Weight: 75 kg
 - Shelves supplied: 2
 - CO2 range/accuracy: 0 to 20/ $\pm 0.1\%$
 - Internal H x W x D: 480 x 560 x 400mm

Didac Scientific Pvt. Ltd,
Martinfield Business Centre, 108 Martinfield, Welwyn Garden City , United Kingdom
Direct Contact Details  sales@didacscientific.co.uk
 www.didacscientific.co.uk/