



Product Code . DS-DASGAGM-10469

Liquid Color Comparator

Description

Liquid Color Comparator

Description:-

The comparator consists of two wheels with nine color filters, ranging from water white through deep amber, imbedded in each wheel. The identifier of the color filter becomes the description for the sample. Simply place a tube of the sample liquid into a slot between the two wheels and rotate the wheels until the color filters are closest in color to the sample.

Didac Scientific Pvt. Ltd,
Martinfield Business Centre, 108 Martinfield, Welwyn Garden City , United Kingdom
Direct Contact Details ✉ sales@didacscientific.co.uk
💻 www.didacscientific.co.uk/