



**Product Code . DS-DASGAGM-10480**

## Advanced Digital Force Gauges

### Description

---



#### Advanced Digital Force Gauges

##### Description:-

These digital force gauges satisfy many force measurement applications in an accurate, sophisticated, and easy-to-use package. The lightweight ergonomic, reversible aluminum design allows for handheld use, or test stand mounting, providing additional flexibility based on your testing requirements. The intuitive menu interface allows for simple navigation and the large backlit display makes access quick and easy to many features and configurable parameters.

- o High speed 3,000 Hz peak sampling rate
- o Live load bar graph with audio/visual overload alarms
- o Backlit graphic display with large digits
- o 50-point data memory with statistics and outputs
- o Peak readings always shown on display
- o Ergonomic, reversible housing for handheld or test stand use
- o Programmable go/no-go limits, with audio/visual alarms and outputs.

---

Direct Contact Details  [sales@didacscientific.co.uk](mailto:sales@didacscientific.co.uk)  
 [www.didacscientific.co.uk/](http://www.didacscientific.co.uk/)