



Product Code . DS-DASGAGM-10483

Professional Digital Force Gauges

Description

Professional Digital Force Gauges



Description:-

Professional digital force gauges are designed with advanced features for your most demanding tension and compression measurement applications, up to 500 lb (2500 N). Gauges feature an industry-leading sampling rate of 14,000 Hz with $\pm 0.1\%$ accuracy, producing quick and reliable results for rapid testing applications.

- o High-speed 14,000 Hz sampling rate
- o Password protection
- o USB, RS-232, Mitutoyo, and analog outputs
- o Sample break detection – automatically zeroes, saves, outputs, and stops a Mark-10 test stand upon sample break
- o Continuous data capture of up to 5,000 data points, at up to 14,000 Hz, downloadable to a PC
- o 1st / 2nd peak detection
- o Five units of measurement, plus COF and user-defined units
- o Programmable footswitch command string

o Averaging, external trigger, real-time, and peak modes

o Programmable set points, with audiovisual indicators and outputs.

Didac Scientific Pvt. Ltd,
Martinfield Business Centre, 108 Martinfield, Welwyn Garden City , United Kingdom
Direct Contact Details  sales@didacscientific.co.uk
 www.didacscientific.co.uk/