



Product Code . DS-DASGAGM-10489

Precision Dynamometer

Description

Precision Dynamometer

Description:-

The Precision Dynamometer is a hand held force gauge for measuring small forces. Small gauges supplied with flat tip and large gauges supplied with round tip. Tolerance ± 0.01 .

Features:-

- o Extreme sensitivity
- o Both clockwise and counter clockwise readings
- o All gauges equipped with maximum reading pointer

Didac Scientific Pvt. Ltd,
Martinfield Business Centre, 108 Martinfield, Welwyn Garden City , United Kingdom
Direct Contact Details ✉ sales@didacscientific.co.uk
💻 www.didacscientific.co.uk/