



Product Code . DS-DASGAGM-10496

Cap Torque Tester

Description

Cap Torque Tester

Features:-

- o Programmable high/low setpoints with LED indicators for quick GO/NO GO testing.
- o Accuracy: $\pm 0.5\%$ F.S. ± 1 LSD
- o Clockwise and counterclockwise operation
- o Clamps containers 1/2" to 8" in diameter
- o Real time and peak mode
- o Tare function
- o AC adapter/charger and carrying case
- o 1000 Data memory for recall and download

Didac Scientific Pvt. Ltd,
Martinfield Business Centre, 108 Martinfield, Welwyn Garden City , United Kingdom
Direct Contact Details ✉ sales@didacscientific.co.uk
💻 www.didacscientific.co.uk/