



**Product Code . DS-DASGAGM-10497**

## Digital Cap Torque Testers

### Description

---



#### Digital Cap Torque Testers

##### Description:-

The digital cap torque testers present a simple, accurate, and economical solution for manual cap torque testing. Testers are ideal for production or lab use for closure manufacturers, bottlers, food and beverage companies, and others. Units capture peak readings for application and removal torques.

- o Mechanical overload protection
  - o Peak memory for application and removal torque readings
  - o RS-232 and Mitutoyo outputs
  - o Accepts samples up to 7.5" (189.4 mm) dia.
  - o Battery or AC operation
  - o Selectable units of measurement (lbin, kgmm, Ncm)
  - o Solid aluminum housing
  - o Moveable sample gripping posts
  - o Accuracy of  $\pm 0.5\%$  of full scale  $\pm 1$  digit
-

---

Didac Scientific Pvt. Ltd,  
Martinfield Business Centre, 108 Martinfield, Welwyn Garden City , United Kingdom  
Direct Contact Details  [sales@didacscientific.co.uk](mailto:sales@didacscientific.co.uk)  
 [www.didacscientific.co.uk/](http://www.didacscientific.co.uk/)