



Product Code . DS-HT-10508

Standard Dial Durometers

Description

Standard Dial Durometers

Description:-

- Standard Dial Durometers feature a large, easy-to-read dial, excellent 1/2 point accuracy, and full 360° sweep.
- Scale is divided into 1-point increments from 0 to 100.
- Measures 2-1/2"W x 6-1/8"H.

Didac Scientific Pvt. Ltd,
Martinfield Business Centre, 108 Martinfield, Welwyn Garden City , United Kingdom
Direct Contact Details ✉ sales@didacscientific.co.uk
💻 www.didacscientific.co.uk/