



Product Code . DS-HT-10509

Test Block Kit

Description

Test Block Kit

Description:-

Test Block Kit for Type A Durometers, 7 Blocks Covering Range 30 to 90. The kit acts as a quick reference check for the status of your durometer and lets you develop the proper feel for use of the gauge. Maintain accurate durometer readings with this kit. Includes seven color-coded rubber test blocks of various type A hardness values ranging from 30 to 90 durometer.

Didac Scientific Pvt. Ltd,
Martinfield Business Centre, 108 Martinfield, Welwyn Garden City , United Kingdom
Direct Contact Details ✉ sales@didacscientific.co.uk
💻 www.didacscientific.co.uk/