



**Product Code . DS-HT-10513**

## Hardness Tester, Ergo

### Description

---

#### Hardness Tester, Ergo

#### Description:-

These instruments measure the indentation hardness of rubber, elastomers, and other rubberlike substances such as neoprene, silicone, vinyl, nitrile, and butyl. Both Durometers come complete with carrying case and test block. The durometers come standard with an ancillary hand for indicating peak readings. They can also be used for soft plastics, felt, leather, and similar materials.

---

Didac Scientific Pvt. Ltd,  
Martinfield Business Centre, 108 Martinfield, Welwyn Garden City , United Kingdom  
Direct Contact Details ✉ [sales@didacscientific.co.uk](mailto:sales@didacscientific.co.uk)  
🌐 [www.didacscientific.co.uk/](http://www.didacscientific.co.uk/)