



Product Code . DS-HT-10515

Impact Testers

Description

Impact Testers

Description:-

- These instruments evaluate impact resistance and determine the exact point of failure or establish pass/fail specifications.
- **Heavy-Duty Tester** is widely used to establish quality control standards for resistance to impact surface damage and penetration of many materials including plastics, resins, fiber glass, sheet metals, and plywood.
- **Light-Duty Tester** is ideal for materials that can be damaged or penetrated by small impact forces, such as products mildly abused in households, offices, or labs through years of normal use.
- **Extra Heavy-Duty Tester** evaluates the impact resistance of rigid sheets of PVC and other materials exceeding the 160 in-lb limit of the heavy-duty tester. Redesigned arm on the base allows high impact to be applied every time.

Didac Scientific Pvt. Ltd,
Martinfield Business Centre, 108 Martinfield, Welwyn Garden City , United Kingdom
Direct Contact Details ✉ sales@didacscientific.co.uk
💻 www.didacscientific.co.uk/