



Product Code . DS-MPA-10520

Beginning/Ending Recording Digital Melting Point Apparatuses

Description

Beginning/Ending Recording Digital Melting Point Apparatuses

Description:-

Programmable Digital Melting Point Apparatus with Beginning/Ending recording for pharmacopeia requirements offers selectable ramp rates between 0.2°C/minute to 10°C/minute with a degree setting of 1°C/minute. Store up to two sets of three temperatures in temporary memory, batch memory permanently stores up to 500 sets of results.

- Digital microprocessor with $\pm 0.1^{\circ}\text{C}$ resolution provides fast warm-up
- Internal batch memory capacity of 500 melts
- Audible beeps indicate when temperature is stable and ready for sample
- Save power with ECO mode
- Adjustable arm with 8x viewing magnifier folds flat for compact storage
- Conveniently located keypad ensures temperatures can be recorded without looking away from the sample.

