



Product Code . DS-MPA-10521

Digital Melting Point Apparatus DID101

Description

Digital Melting Point Apparatus DID101

Description:-

Designed for maximum safety and ease of use. Select, measure, and display the temperature digitally, ensuring accuracy and avoiding the need for a glass thermometer. Extend the back feet for the optimal viewing angle. Set a plateau temperature for initial warm-up — unattended warm-up time frees up the user until the plateau temperature is reached. Test two samples simultaneously, viewing the melt via a magnifying lens with built-in illumination. Heating rate continues at 2°C after the warm-up, slowly reaching the melt temperature for greater accuracy.

- Select, measure, and display the temperature digitally, ensuring accuracy and repeatability
- Save time by setting a plateau temperature for initial warm-up
- Be more productive—test multiple samples simultaneously
- Easily view the melt via magnifying lens with built-in illumination.

Didac Scientific Pvt. Ltd,
Martinfield Business Centre, 108 Martinfield, Welwyn Garden City , United Kingdom
Direct Contact Details ✉ sales@didacscientific.co.uk
🌐 www.didacscientific.co.uk/