



**Product Code . DS-MPA-10522**

## **Standard Digital Melting Point Apparatus**

### **Description**

---

#### **Standard Digital Melting Point Apparatus**



##### **Description:-**

Digital Melting Point Apparatus is ideal for basic melting point testing of up to three samples. Use the 10°C/minute search ramp to quickly analyze sample if the melt temperature is unknown. A fixed 1°C/minute ramp heats sample. Store up to four temperatures in temporary memory.

##### **Key Features:-**

- Digital microprocessor with  $\pm 0.1^{\circ}\text{C}$  resolution provides fast warm-up
  - Save power with ECO mode
  - Stores up to four temperatures in temporary memory
  - Adjustable arm with 8x viewing magnifier folds flat for compact storage
  - Audible beeps indicate when temperature is stable and ready for sample
  - Conveniently located keypad ensures temperatures can be recorded without looking away from the sample.
-

---

Didac Scientific Pvt. Ltd,  
Martinfield Business Centre, 108 Martinfield, Welwyn Garden City , United Kingdom  
Direct Contact Details  [sales@didacscientific.co.uk](mailto:sales@didacscientific.co.uk)  
 [www.didacscientific.co.uk/](http://www.didacscientific.co.uk/)