



**Product Code . DS-MPA-10523**

## Mel-Temp Melting Point Apparatus

### Description

---

#### Mel-Temp Melting Point Apparatus

##### Description:-

The Melting Point temperature of up to three samples using the Mel-Temp apparatus. Simply key in your set point temperature and once the oven temperature reaches the set temperature, audible beeps indicate oven temperature stability. Position the adjustable arm where you need it and rotate the viewing head to eliminate glare on the lens. The temperature slowly ramps up from the set point and the melt is observed through the 8x magnifying lens, melt temperature is displayed on the LCD. Up to four melt temperatures can be stored in temporary memory.

Power (VAC): 100 to 240

Power (Hz): 50/60

---

Didac Scientific Pvt. Ltd,  
Martinfield Business Centre, 108 Martinfield, Welwyn Garden City , United Kingdom  
Direct Contact Details ✉ [sales@didacscientific.co.uk](mailto:sales@didacscientific.co.uk)  
💻 [www.didacscientific.co.uk/](http://www.didacscientific.co.uk/)