



Product Code . DS-MPA-10524

Programmable Digital Melting Point Apparatus

Description



Programmable Digital Melting Point Apparatus

Description:-

Programmable ramp rate digital melting point apparatus features five preset temperature ramps, or you can program for any rate from 0.2°C/minute to 10°C/minute in 0.1°C increments. Store up to four temperatures in temporary memory, batch memory permanently stores up to 1000 sets of results with date, time, and batch number. Power 100-240 VAC, 50-60 Hz.

Key Features:-

- Digital microprocessor with $\pm 0.1^{\circ}\text{C}$ resolution provides fast warm-up
 - Conveniently located keypad ensures temperatures can be recorded without looking away from the sample
 - Internal batch memory capacity of 1000 melts
 - Save power with ECO mode
 - Adjustable arm with 8x viewing magnifier folds flat for compact storage
 - Audible beeps indicate when temperature is stable and ready for sample.
-

Didac Scientific Pvt. Ltd,
Martinfield Business Centre, 108 Martinfield, Welwyn Garden City , United Kingdom
Direct Contact Details  sales@didacscientific.co.uk
 www.didacscientific.co.uk/