



Product Code . DS-MPA-10525

Advanced Digital Melting Point Apparatus

Description

Advanced Digital Melting Point Apparatus

Description:-

Advanced Digital Melting Point Apparatus features an innovative “head-up” display which projects a floating red image of the block temperature in the same viewing field as the magnified tubes. In addition, the viewing tube can be adjusted up or down for a comfortable viewing angle and then folds flat for easy storage. The advanced unit lets you test three samples simultaneously and features a maximum temperature of 400°C, a variable ramp, and three programmable methods for repeated use. This integrated display eliminates the need for you to constantly switch your gaze between the tubes and the temperature shown on the control panel. Large backlit LCD shows programming information, stored results, current temperatures, and time. Up to seven date-and-time stamped results per tube can be stored in memory.

Didac Scientific Pvt. Ltd,
Martinfield Business Centre, 108 Martinfield, Welwyn Garden City , United Kingdom
Direct Contact Details ✉ sales@didacscientific.co.uk
💻 www.didacscientific.co.uk/