

Product Code . DS-MPA-10528

Analog Melting Point Apparatuses DID102



Description

Analog Melting Point Apparatuses DID102

Description:-

The analog melting point apparatus features a blue-spirit-filled thermometer for mercury-free temperature measurements. The ramp rate is 1 to 20°C per minute until melt occurs. Test and view up to three samples at once for greater efficiency. The unit is simple to operate with instructions printed clearly on the front panel — ideal for educational settings.

- Easily view the melt via magnifying lens with built-in illumination.
- Be more productive—test multiple samples simultaneously.

Didac Scientific Pvt. Ltd,
Martinfield Business Centre, 108 Martinfield, Welwyn Garden City , United Kingdom
Direct Contact Details ✉ sales@didacscientific.co.uk
🌐 www.didacscientific.co.uk/