



**Product Code . DS-MPA-10530**

## **Fluorescent Indicator Adsorption Apparatus**

### **Description**

---

#### **Fluorescent Indicator Adsorption Apparatus**

##### **Description:-**

Optimize performance of Fluorescent Indicator Adsorption Apparatus Compatible for use with the Fluorescent Indicator Absorption Apparatus, sold separately. After the sample is introduced and adsorbed into the column, fluorescent dye indicators mark the boundaries of aromatics, olefins, and saturates, which are visible with a UV lamp. An integrated electric vibration system facilitates the dry gel packing of the columns.

---

Didac Scientific Pvt. Ltd,  
Martinfield Business Centre, 108 Martinfield, Welwyn Garden City , United Kingdom  
Direct Contact Details ✉ [sales@didacscientific.co.uk](mailto:sales@didacscientific.co.uk)  
💻 [www.didacscientific.co.uk/](http://www.didacscientific.co.uk/)