



**Product Code . DS-MPA-10532**

## **Ramsbottom Carbon Residue Furnace**

### **Description**



---

#### **Ramsbottom Carbon Residue Furnace**

##### **Description:-**

Heavily insulated stainless steel cabinet with three-layer refractory top provides excellent heat retention. The microprocessor temperature control keeps 2°F (1°C) stability. Determine the carbon residue of your sample which provides an indication of the deposit-forming tendencies of fuels. Simply insert up to five samples, the furnace quickly reaches the test temperature of 1022°F (550°C). Apparatus features digital temperature setting, over temperature control, ceramic band heater with cast iron block, and heavily insulated stainless steel cabinet.

---

Didac Scientific Pvt. Ltd,  
Martinfield Business Centre, 108 Martinfield, Welwyn Garden City , United Kingdom  
Direct Contact Details  [sales@didacscientific.co.uk](mailto:sales@didacscientific.co.uk)  
 [www.didacscientific.co.uk/](http://www.didacscientific.co.uk/)