



**Product Code . DS-MPA-10533**

## **Salt-In-Crude Analyzer**

### **Description**

---

#### **Salt-In-Crude Analyzer**

##### **Description:-**

Quickly analyze oil samples with this salt-in-crude analyzer and make the necessary operation changes to prevent corrosion of equipment and pipe lines due to excessive salt content. The self-calibration feature allows the operator to adjust for any drift without re-entering standard temperature curves. Analyzer measures conductance, and temperature (°C or °F) of crude oil samples, and pH of aqueous samples. Automatically calculates salt content from measured conductance and temperature values. Measure conductivity over four ranges 0-2, 2-20, 20-200, and 200-1500  $\mu$ S.

---

Didac Scientific Pvt. Ltd,  
Martinfield Business Centre, 108 Martinfield, Welwyn Garden City , United Kingdom  
Direct Contact Details ✉ [sales@didacscientific.co.uk](mailto:sales@didacscientific.co.uk)  
🌐 [www.didacscientific.co.uk/](http://www.didacscientific.co.uk/)