



Product Code . DS-MPA-10535

Copper Corrosion Test Bath System

Description

Copper Corrosion Test Bath System

Description:-

Determine the corrosiveness to copper of petroleum-based fuels and bio diesel by immersing a copper strip in the test sample at an elevated temperature. Built-in support rack holds up to 16 test tubes. Over temperature control interrupts power to the heater if programmed temperature is reached. The constant-temperature bath features a digital temperature controller and a 750-watt heater, maintains a temperature stability of $\pm 2^{\circ}\text{F}$ (1°C).

Didac Scientific Pvt. Ltd,
Martinfield Business Centre, 108 Martinfield, Welwyn Garden City , United Kingdom
Direct Contact Details ✉ sales@didacscientific.co.uk
💻 www.didacscientific.co.uk/