



Product Code . DS-MPA-10537

Existent Gum Test Apparatus

Description

Existent Gum Test Apparatus

Description:-

Steam Super heater provides dry super heated steam for evaporation baths not equipped with a built-in super heater. Super heating chamber and condensate trap are constructed entirely of stainless steel. Use together with an outside steam source for steam-jet method testing of aircraft turbine fuels. Solid state temperature controller adjusts between 0-550°F. Steel exterior has a chemical resistant polyurethane enamel finish. Equipped with steam inlet and outlet connections and condensate drain valve.

Didac Scientific Pvt. Ltd,
Martinfield Business Centre, 108 Martinfield, Welwyn Garden City , United Kingdom
Direct Contact Details ✉ sales@didacscientific.co.uk
🌐 www.didacscientific.co.uk/