



**Product Code . DS-MPA-10538**

## **Pensky-Martens Gas-Heated Closed-Cup Flash Tester**

### **Description**

---

#### **Pensky-Martens Gas-Heated Closed-Cup Flash Tester**

##### **Description:-**

Pensky-Martens Closed-Cup Flash Testers feature a cover that slides shutter open to apply test flame. Electrical models are equipped with a 750 watt nickel-chromium heater. Gas model has a built-in nickel-plated brass natural gas burner and liquid propane burner mounted on a cast iron base. Select from electrical or gas heated models. Testers include a thermometer holder, brass test cup with handle, pilot flame, test flame reference bead, built-in stirrer, plated brass thermometer ferrule, and cooling vents.

---

Didac Scientific Pvt. Ltd,  
Martinfield Business Centre, 108 Martinfield, Welwyn Garden City , United Kingdom  
Direct Contact Details ✉ [sales@didacscientific.co.uk](mailto:sales@didacscientific.co.uk)  
💻 [www.didacscientific.co.uk/](http://www.didacscientific.co.uk/)