

Product Code . DS-MPA-10541

Automatic Polarimeter



Description

Automatic Polarimeter

Description:-

Take automatic measurements of angle of rotation and International Sugar Scale in samples such as sugars, food additives, pharmaceuticals, and chemicals. Adjust operational parameters such as scale, tube length, and temperature through the menu-driven touch screen. Download your results to a computer or printer through the RS-232 interface. Perform an instrument blank with an empty tube — sample readings are automatically adjusted based on the blank. Select between English and German language for the display. Power: 100 to 250 VAC, 50/60 Hz.

Didac Scientific Pvt. Ltd,
Martinfield Business Centre, 108 Martinfield, Welwyn Garden City , United Kingdom
Direct Contact Details ✉ sales@didacscientific.co.uk
🌐 www.didacscientific.co.uk/