

**Product Code . DS-MPA-10543**



## Two-Gauge Manifold for Manual Hydraulic Press Test

### Description

---

#### Two-Gauge Manifold for Manual Hydraulic Press Test

##### Description:-

Hand-operated, two-gauge manifold includes a shut-off valve for the low-pressure gauge to protect it during high-pressure applications. Two-Gauge Manifold provides both low- and high-pressure gauge readings for greater precision. Unit includes a connection for two gauges.

---

Didac Scientific Pvt. Ltd,  
Martinfield Business Centre, 108 Martinfield, Welwyn Garden City , United Kingdom  
Direct Contact Details ✉ [sales@didacscientific.co.uk](mailto:sales@didacscientific.co.uk)  
💻 [www.didacscientific.co.uk/](http://www.didacscientific.co.uk/)