



**Product Code . DS-MPA-10544**

## **Hydraulic Press Test Systems, Automatic**

### **Description**

---

#### **Hydraulic Press Test Systems, Automatic**

##### **Description:-**

Automatic Hydraulic Presses feature a touch-screen interface allowing you to select auto or manual mode, and units in lb, US tons, metric tons, or kg. Presses feature audible/visual alarms and proximity switch for adjustable slowdown position. Requires two-hand operation for maximum safety. Scroll through menu to enter parameters for pressure, pump speed, and dwell time. Safety features include a 10-gauge steel shield with see-through, shatterproof polycarbonate access door. These heavy-duty systems handle a wide range of laboratory quality and performance applications such as compression, flexural, or shear strength, flow and color dispersion, and material research using ceramics, composites, construction materials, resins, and more. The base has holes for permanent bench mounting and for optional permanent bench mounting. Compact design efficiently uses your laboratory bench space. Cast-iron base and bolsters with two threaded columns provide easy daylight adjustment.

---

Didac Scientific Pvt. Ltd,  
Martinfield Business Centre, 108 Martinfield, Welwyn Garden City , United Kingdom  
Direct Contact Details ✉ [sales@didacscientific.co.uk](mailto:sales@didacscientific.co.uk)  
💻 [www.didacscientific.co.uk/](http://www.didacscientific.co.uk/)