



Product Code . DS-MPA-10545

Hydraulic Press Test System, Manual

Description

Hydraulic Press Test System, Manual

Description:-

These heavy-duty systems handle a wide range of laboratory quality and performance applications such as compression, flexural, or shear strength, flow and color dispersion, and material research using ceramics, composites, construction materials, resins, and more. Pressure gauge indicates both applied force and maximum pressure. Handle has cushioned comfort grip. Safety features include a 10-gauge steel shield with see-through, shatterproof poly carbonate access door. The base has holes for permanent bench mounting and for optional permanent bench mounting.

Didac Scientific Pvt. Ltd,
Martinfield Business Centre, 108 Martinfield, Welwyn Garden City , United Kingdom
Direct Contact Details ✉ sales@didacscientific.co.uk
💻 www.didacscientific.co.uk/