

**Product Code . DS-SAS-10549**

## UV/Visible Spectrophotometers



### Description

---

#### UV/Visible Spectrophotometers

##### Description:-

UV/Visible spectrophotometers have a low profile with a large LED display which can be read from any angle. The 4 nm bandwidth of these single-beam instruments produces precise response to transmittance changes for UV/visible spectrum methods. Four-position cuvette holders save time and improve efficiency. Use the USB ports for connecting directly to a printer or to a PC. Peltier controller regulates temperature between 15 and 40°C.

- Software-driven DNA/RNA/protein/kinetics analysis.
- Store up to 200 test methods for your demanding lab environments.

---

Didac Scientific Pvt. Ltd,  
Martinfield Business Centre, 108 Martinfield, Welwyn Garden City , United Kingdom  
Direct Contact Details ✉ [sales@didacscientific.co.uk](mailto:sales@didacscientific.co.uk)  
💻 [www.didacscientific.co.uk/](http://www.didacscientific.co.uk/)