

Product Code . DS-SAS-10554

Benchtop Visible Spectrophotometers DID103



Description

Benchtop Visible Spectrophotometers DID103

Description:-

This easy-to-use spectrophotometer is ideal for the busy technician or student. Wavelength calibration is automatic at start-up and includes an auto-zero function—a significant time saver. The included 63-Zero PC software allows you to take and control measurements remotely, and export data and results to Excel, Word, Power Point, and other software applications for storage, reporting, or presentation. Transmittance, absorbance, or concentration with factor value, as well as wavelength are simultaneously. This spectrophotometer is fully GLP compliant. Each sample reading is given a time and date stamp and easily connects to a printer—then simply push the print key to send data. The digital, dual-line LCD shows error messages, prompts, mode indication, and concentration units.

- Dual-line LCD shows error messages, prompts, mode indication, and concentration units.
- Recalls previous settings at power up.
- Analog and RS-232 outputs let you interface with printers and PCs.

Didac Scientific Pvt. Ltd,
Martinfield Business Centre, 108 Martinfield, Welwyn Garden City , United Kingdom
Direct Contact Details ✉ sales@didacscientific.co.uk
💻 www.didacscientific.co.uk/