



Product Code . DS-SAS-10560

Nano Photometer UV/Vis Spectrophotometers

Description

Nano Photometer UV/Vis Spectrophotometers

Description:-

For those in need of flexibility, this family of UV/visible spectrophotometers offers reliable protein and DNA measurements for mobile and laboratory applications. Simplified optics with two precisely defined pathlengths provide two fixed and sealed anchor points, prohibiting drift and eliminating the need for recalibration. Requiring only a minute sample size of 0.3 μL and with a limit of detection down to 0.1 ng/ μL (ds DNA) or 0.003 mg/ μL (BSA), this spectrophotometer series accommodates any application. Equipped with Wi-Fi, USB A/B, HDMI and LAN interface connections for flexible instrument control and data saving options. Front USB flash drive port offers quick and easy data transfer.

- Flexible control and data collection via computer, Mac, smartphone, or tablet.
- Choose a Touch model for a built-in touch screen tablet.

Didac Scientific Pvt. Ltd,
Martinfield Business Centre, 108 Martinfield, Welwyn Garden City , United Kingdom
Direct Contact Details ✉ sales@didacscientific.co.uk
💻 www.didacscientific.co.uk/