

Product Code . DS-SAS-10561

UV-Vis Spectrophotometer DID104



Description

UV-Vis Spectrophotometer DID104

Description:-

This spectrophotometer revolutionizes instrument calibration with optional RFID technology. This results in a more efficient, accurate instrument. On top of that, the unit has 250 preprogrammed methods ready to go. The RFID model will recognize a TNTplus test when placed near the scanner, automatically comparing your individual calibration information to the factory calibration curve. This makes running any common test parameter easy and painless. The units have flexible connectivity and come with one ethernet port and three USB ports.

- Operates in the UV and visible spectrum
- Features over 250 pre-programmed procedures for water analysis
- The instrument takes 10 readings of one sample to reduce error and remove outliers.

Didac Scientific Pvt. Ltd,
Martinfield Business Centre, 108 Martinfield, Welwyn Garden City , United Kingdom
Direct Contact Details ✉ sales@didacscientific.co.uk
💻 www.didacscientific.co.uk/