

Product Code . DS-SAS-10562

Spectrophotometer



Description

Spectrophotometer

Description:-

These spectrophotometers revolutionize instrument calibration with optional RFID technology. The RFID model will recognize a TNTplus test when placed near the scanner, automatically comparing your individual calibration information to the factory-installed calibration curve. The units have flexible connectivity and come with one ethernet port and three USB ports. This results in a more efficient, accurate instrument. On top of that, the unit has 220 preprogrammed methods ready to go. This makes running any common test parameter easy and painless.

- Feature over 220 preprogrammed procedures for water analysis
- Options available for EPA compliant test parameters
- The instrument takes 10 readings of one sample to reduce error and remove outliers.

Didac Scientific Pvt. Ltd,
Martinfield Business Centre, 108 Martinfield, Welwyn Garden City , United Kingdom
Direct Contact Details ✉ sales@didacscientific.co.uk
🌐 www.didacscientific.co.uk/